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6-Axis Collaborative Robot

Warning

Before using myCobot 320 Pi, please read all instructions and cautionary markings in this manual

1. Do not expose the product to rain or moisture to reduce fire or shock hazard.
2. Do not place the product in or near fire.
3. Do not leave the product in a car in hot or humid weather.
4. Do not disassemble, crush or pierce the product.
5. Do not expose the product to excessive shock such as dropping from a high place.
6. Do not expose the product to high temperatures above 60 °C (140 °F).

Attention


Regarding the operation and secondary development of myCobot Pro Pi, please read and download **Gitbook** before using it.

Official Website: <https://www.elephantrobotics.com/en/mycobot-pro-pi-en/>

Model: myCobot-320-Pi

myCobot 320 Pi


Highly Cost-Effective

6-Axis Collaborative Robot

myCobot 320 Pi is a product belonging to the "myCobot 320 series". It uses a Raspberry Pi micro-processor and embedded RoboFlow visual programming software. It is a machine "assistant" developed by Elephant Robotics for makers and researchers.

The myCobot 320 Pi has a weight of 3kg, a load of 1kg, and a working radius of 320mm. The volume is relatively small but powerful. It is easy to operate, can work in collaboration with people, and work safely. It has three advantages: ease of use, safety and economy, and it is a cost-effective choice.

Embedded Raspberry Pi ecology, unlimited development possibilities

- Raspberry Pi 4B, 1.5GHz 4-core microprocessor, running on Debian/Ubuntu platform.
- Support 4 USB, 2 HDMI, standardized GPIO interface, TF card pluggable.

Easy to operate, open source and easy to use

- Drag to teach, users can quickly get started.
- Support development systems such as ROS/moveIt and RoboFlow operating software independently developed by Elephant Robotics.

High configuration, powerful performance

- The use of brushless DC servos can achieve a repeat positioning accuracy of $\pm 0.5\text{mm}$.
- The base and the end are equipped with installation interfaces, which are suitable for the development of various peripheral products and equipment.

Economical and applicable, super cost-effective

- The standard 8-hour work system can replace repetitive and standard work.
- Cost effective for reducing costs and increasing efficiency for scientific research that requires high performance and low cost.

Integrated design, safe collaborative operation

- The exquisite structural design makes it possible to make full use of the space and perfectly blend into the actual environment.
- Equipped with the anti-collision detection function based on the accurate dynamic model, so that it can work safely with people.


Unique industrial design


- Integrated design, the overall structure is compact, and the net weight is 3kg.
- Modular design, few spare parts, low maintenance cost, quick disassembly and replacement.

myCobot 320 - Design Prototype - Elephant Robot®C Series Robot

The design prototype of myCobot 320 is from All-in-one Robot launched by Elephant Robot in China in 2018. As the first integrated collaborative robot in China, it has won the 2019 CAIMRS Industrial Robot Innovation Award and 2019 High-tech Robot Annual "Innovation Technology Award", and has been also sold to more than 30 countries at home and abroad, receiving unanimous praise and recognition from the factories of the world's top 500 enterprises.

myCobot 320 Pi - Size and Working Range Diagram

Parameter

| Robot arm Parameter | |
|----------------------------|--------------------------|
| Model | myCobotPro 320-Pi |
| Degree of Freedom | 6 |
| Payload | 1kg |
| Working range | 320mm |
| Repeatability | $\pm 0.5\text{mm}$ |
| Weight | 3kg |
| Power Input | 24V, 5A |
| Working Temp. | 0-50° |
| Communication | Type-C |
| Software Parameter | |
| Programming Platform | Debian/Ubuntu |
| ROS/Python | Built-in |
| Blockly Visual Programming | Built-in |
| Camera | Available for any camera |
| ModelBus RTU | Ethernet |

| Electric parameter | |
|------------------------------|-------------------------|
| SOC | Broadcom BCM2711 |
| CPU | 64-bit 1.5GHz quad-core |
| Bluetooth/wireless | Yes |
| USB | USB3.0 x2; USB2.0 x2 |
| Display | No |
| HDMI port | microHDMI x2 |
| Customizable button | No |
| IO Port | 8 |
| Emergency Stop Button | Yes |
| Teach pendant/Keyboard/Mouse | Optional |

Control Board Pin Map

Atom diagram

Atom pin diagram

myStudio

myStudio is a one-stop platform for robots

myStudio integrates myCobot's software and various materials. The main functions of myStudio are: 1) Update the firmware; 2) Provide video tutorials on how to use the robot; 3) Provide maintenance and repair information (such as video tutorials, Q&A, etc.).

Please download the latest version of myStudio to use.


The download link is as follows:

Official website: <https://www.elephantrobotics.com/myCobot/>

Github: <https://github.com/elephantrobotics/MyStudio/>

Burn Table

Development environments that support the secondary development of myCobot 320 are: Myblockly, RoboFlow, Arduino, ROS, python, etc.

| Development Environment | Library on PC | Atom Firmware |
|---|--|---------------|
| Default Program | N/a | atomMain |
| Myblockly | Myblockly | atomMain |
| RoboFlow Industrial Programming Software | RoboFlow | atomMain |
| Arduino Maker! | Arduino IDE + M5Stack Lib + MycobotBasic Lib | atomMain |
| API on Desktop | Python/ C++ | atomMain |
| ROS Development | ROS | atomMain |
| USB/TxRx0(G1/G3) | Read Protocol | atomMain |
| BlueTooth | Read Protocol | atomMain |
| phoneApp | Mobile phone Android/iPhone | atomMain |


myCobot 320 Accessory

Elephant Robotics are targeted at robotic collaboration applications, making “my-series” product line. For new information about the accessories, Follow us on Shopify and Twitter.

Shopify: <https://shop.elephantrobotics.com/>

Twitter: @cobotMy

Facebook: myCoboy

- 1 Electric parallel gripper
- 2 Air parallel gripper
- 3 Adaptive gripper
- 4 Vacuum suction cups

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Purchaser _____ Order No. _____ Phone _____

Address _____ Logistics Receipt Date _____

Product problem description(Required):

If you need to apply for warranty service, please contact our customer service to confirm the detailed information. After confirmation, please fill in the card and send it back together with the product and the attached invoice. **Note: Our company reserves the right to explain and modify the warranty card of this product within the scope of the law.**

- Return service is limited to goods not opened within 7 days after the receipt date of logistics of the products. The freight or other risks incurred in return shall be borne by the customer.
- Customers should provide the purchasing invoice and warranty card as the warranty certification when a warranty is being asked.
- Elephant Robotics will be responsible for the hardware faults of products caused by the normal using during the warranty period.
- The warranty period starts from the date of purchase or the receipt date of the logistics.
- The faulty parts from the products will be owned by Elephant Robotics, and the appropriate cost will be charged if necessary.

If you need to apply for warranty service, please contact our customer service first to confirm the detailed information.

Sever motor

Warranty Period Warranty Services

| | |
|------------|--|
| ≤1 months | Elephant Robotics offers a free new servo motor and bear the freight. |
| 1-3 months | Elephant Robotics offers a free new servo motor, customs shall bear the freight. |
| ≥3 months | Customers need to buy it themselves. |

Electrical Parts

| | |
|------------|--|
| ≤3 months | Customers need to send it back after disassembly, Elephant Robotics shall send a new one for free and bear the freight out and home. |
| 3-6 months | Customers need to send it back after disassembly and bear the freight out and home, Elephant Robotics shall send a new one for free. |
| ≥6 months | Customers need to buy it themselves. |

Structure Parts, including Shell Parts

| | |
|---------|--|
| ≤1 year | Elephant Robotics offers free new components once, customs shall bear the freight. |
| ≥1 year | Customers need to buy it themselves. |

During the warranty period of the delivered product, the company only repairs the malfunctions that occur during normal use of the robot for free. However, in the following cases, the customer will be charged for repairs (even during the warranty period):

- Damage or malfunction caused by incorrect use and improper use different from the contents of the manual.
- Failure caused by unauthorized disassembly by the customer.
- Damage caused by improper adjustment or unauthorized repairs.
- Damage caused by natural disasters such as earthquakes and floods.

Therefore, please strictly follow the instructions in this manual and related manual to operate the robot.

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网址：www.elephantrobotics.com

Shenzhen Elephant Robotics Technology Co., Ltd

Address: B7, Yungu Innovative Industrial Park 2, Nanshan, Shenzhen, China

Email: support@elephantrobotics.com

Phone: +86(0755)-8696-8565 (working day 9:30-18:30)

Website: www.elephantrobotics.com